



## CPR Feedback Change Notice



### The service is **phasing in** a new CPR feedback technology

You may use 2 different technologies of CPR feedback sensors until the rollout has been completed. Below are the differences between the cables and CPR feedback sensors for the **old** and **new** technology setup. They are **not interchangeable**.

#### Phasing Out Old Technology corPatch CPR Feedback

Identified by **BLUE** connectors



#### Phasing In New Technology Reusable PrimeCPR Feedback

Identified by **GREEN** connectors



#### Ortus Academy

Ortus Academy is a trading name of Ortus Technology Limited.  
Company Registration 07554900